



# Liposuction Cannulas & Accessories



### **Table of Contents**

- A. Liposuction Cannulas
  - a) Infiltration Cannulas
  - b) General Suction Cannulas
  - c) Harvesting Cannulas
  - e) Injecting Cannulas
  - f) Speciality Cannulas
- **B.** Accessories
- C. Our Specialized Sets



### Cannula Hubs/ Connections







Luer lock



Super Luer Lock



Treaded



60cc Tommey hub

### **Machine Connections**





### **Luer Lock Hub/Connector:**

The universal Luer Lock hub enables cannula to attach directly with any Luer Lock Syringe.

### **Super Luer Lock Hub/Connector:**

A modification to general Luer Lock hub which is more powerful, more stable, Leak proof, free from instability and contamination.

### **Threaded Hub/Connector:**

Threaded connector hub style is used mostly for machine assisted Liposuction cannulas, Medicta's unique Light weight threaded cannula can also be used with ordinary syringe. All you need is Medicta's connector from threaded to Luer style.

### **60cc Tommey Hub/Connector:**

An ideal hub style for surgeons why love to handle large fat volumes. This hub style connect cannula to a 60cc Tommey to syringe.

### **Fixed Handle Cannulas:**

The unique design of handle offer powerful grip harvesting large volumes. You will feel amazingly comfortable ever after a long operational time.

The thumb grip provide a unique opportunity to control the thrust.



### INFILTRATION CANNULAS

Infiltration Cannulas are designed for safe and theicfifedelivery of tumescent solution. The spiral micro holes on cannula ensure and even distribution of local anesthetic which results in a longer anesthetic effects this less painful and more pleasent liposuction experience for patient.



### **INFILTRATION CANNULAS**



#### **Tumescent Infiltrator**

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

	60cc Tommey hub
ML-101.01 ML-101.02 ML-101.03 ML-101.04	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm
	Treaded Hub
ML-102.01 ML-102.02 ML-102.03 ML-102.04 ML-102.05 ML-102.06	Ø 2.0 mm X 15 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 30 cm Ø 4.0 mm X 20 cm Ø 4.0 mm X 30 cm
	Super Luer Lock
ML-103.01	Ø 1.0 mm X 5 cm

	Luer lock	
ML-104.01	Ø 1.0 mm X 4 cm	
ML-104.02	Ø 1.0 mm X 5 cm	
ML-104.03	Ø 1.5 mm X 7 cm	
ML-104.04	Ø 1.5 mm X 9 cm	
ML-104.05	Ø 1.5 mm X 10 cm	
ML-104.06	Ø 1.5 mm X 12 cm	
ML-104.07	Ø 2.0 mm X 10 cm	
ML-104.08	Ø 2.0 mm X 12 cm	
ML-104.09	Ø 2.0 mm X 15 cm	
ML-104.10	Ø 2.0 mm X 20 cm	
ML-104.11	Ø 2.5 mm X 15 cm	
ML-104.12	Ø 2.5 mm X 20 cm	
ML-104.13	Ø 3.0 mm X 15 cm	
MI -104 14	Ø 3.0 mm X 20 cm	

Ø 1.0 mm X 7 cm Ø 1.5 mm X 7 cm

Ø 1.5 mm X 9 cm

	Luer lock (Stainless Steel)
ML-105.01	Ø 1.0 mm X 4 cm
ML-105.02	Ø 1.0 mm X 5 cm
ML-105.03	Ø 1.0 mm X 7 cm
ML-105.04	Ø 1.5 mm X 7 cm
ML-105.05	Ø 1.5 mm X 9 cm
ML-105.06	Ø 1.5 mm X 10 cm
ML-105.07	Ø 1.5 mm X 12 cm
ML-105.08	Ø 2.0 mm X 10 cm
ML-105.09	Ø 2.0 mm X 12 cm
ML-105.10	Ø 2.0 mm X 15 cm
ML-105.11	Ø 2.0 mm X 20 cm
ML-105.12	Ø 2.5 mm X 15 cm
ML-105.13	Ø 2.5 mm X 20 cm
ML-105.14	Ø 3.0 mm X 15 cm
ML-105.15	Ø 3.0 mm X 20 cm

	Power Cannula	
ML-106.01	Ø 2.0 mm X 15 cm	
ML-106.02	Ø 2.0 mm X 20 cm	
ML-106.03	Ø 2.0 mm X 25 cm	
ML-106.04	Ø 2.0 mm X 30 cm	
ML-106.05	Ø 2.0 mm X 35 cm	
ML-106.06	Ø 3.0 mm X 15 cm	
ML-106.07	Ø 3.0 mm X 20 cm	
ML-106.08	Ø 3.0 mm X 25 cm	
ML-106.09	Ø 3.0 mm X 30 cm	
ML-106.10	Ø 3.0 mm X 35 cm	
ML-106.11	Ø 4.0 mm X 15 cm	
ML-106.12	Ø 4.0 mm X 20 cm	
ML-106.13	Ø 4.0 mm X 25 cm	
ML-106.14	Ø 4.0 mm X 30 cm	
ML-106.15	Ø 4.0 mm X 35 cm	
ML-106.16 ML-106.17	Ø 5.0 mm X 15 cm Ø 5.0 mm X 20 cm	
ML-106.17	Ø 5.0 mm X 25 cm	
ML-106.19	Ø 5.0 mm X 30 cm	
ML-106.19	Ø 5.0 mm X 35 cm	
ML-106.21	Ø 6.0 mm X 15 cm	
ML-106.21	Ø 6.0 mm X 20 cm	
ML-106.23	Ø 6.0 mm X 25 cm	
ML-106.24	Ø 6.0 mm X 30 cm	
ML-106.25	Ø 6.0 mm X 35 cm	
WIE 100.20	2 0.0 11111 7 00 0111	

ML-103.02 ML-103.03 ML-103.04







Tumescent Bayonet Infiltrator
The bayonet shaped handles allow an easier manipulation

	60cc Tommey hub	
ML-101.11	Ø 2.0 mm X 15 cm	
ML-101.12	Ø 2.5 mm X 20 cm	
ML-101.13	Ø 3.0 mm X 20 cm	
ML-101.14	Ø 3.0 mm X 25 cm	


	Treaded Hub
ML-102.11	Ø 2.0mm X 15 cm
ML-102.12	Ø 3.0mm X 15 cm
ML-102.13	Ø 3.0mm X 20 cm
ML-102.14	Ø 3.0mm X 30 cm
ML-102.15	Ø 4.0mm X 20 cm
WIL-102.10	
	Super Luer Lock
ML-103.11	Ø 1.0 mm X 5 cm
	~
ML-103.12	Ø 1.0 mm X 7 cm
ML-103.13	Ø 1.5 mm X 7 cm
ML-103.14	Ø 1.5 mm X 9 cm

	Luer lock
ML-104.11	Ø 1.0 mm X 4 cm
ML-104.12	Ø 1.0 mm X 5 cm
ML-104.13	Ø 1.0 mm X 7 cm
ML-104.14	Ø 1.5 mm X 7 cm
ML-104.15	Ø 1.5 mm X 9 cm
ML-104.16	Ø 1.5 mm X 10 cm
ML-104.17	Ø 1.5 mm X 12 cm
ML-104.18	Ø 2.0 mm X 10 cm
ML-104.19	Ø 2.0 mm X 12 cm
ML-104.20	Ø 2.0 mm X 15 cm
ML-104.21	Ø 2.0 mm X 20 cm
ML-104.22	Ø 2.5 mm X 15 cm
ML-104.23	Ø 2.5 mm X 20 cm
ML-104.24	Ø 3.0 mm X 15 cm
ML-104.25	Ø 3.0 mm X 20 cm

	Luer lock (Stainless Steel)	
ML-105.21	Ø 1.0 mm X 4 cm	
ML-105.22 ML-105.23	Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm	
ML-105.23	Ø 1.5 mm X 7 cm	
ML-105.25	Ø 1.5 mm X 9 cm	
ML-105.26	Ø 1.5 mm X 10 cm	
ML-105.27	Ø 1.5 mm X 12 cm	
ML-105.28	Ø 2.0 mm X 10 cm	
ML-105.29	Ø 2.0 mm X 12 cm	
ML-105.30	Ø 2.0 mm X 15 cm	
ML-105.31	Ø 2.0 mm X 20 cm	
ML-105.32	Ø 2.5 mm X 15 cm	
ML-105.33 ML-105.34	Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm	
IVII = 1115 34	Power Cannula	
MI 400.04		
ML-106.31	Ø 2 mm X 15 cm Ø 2 mm X 20 cm	
ML-106.32 ML-106.33	Ø 2 mm X 25 cm	
ML-106.34	Ø 2 mm X 30 cm	
ML-106.35	Ø 2 mm X 35 cm	
ML-106.36	Ø 3 mm X 15 cm	
ML-106.37	Ø 3 mm X 20 cm	
ML-106.38	Ø 3 mm X 25 cm	
ML-106.39	Ø 3 mm X 30 cm	
ML-106.40	Ø 3 mm X 35 cm	
ML-106.41	Ø 4 mm X 15 cm	
ML-106.42	Ø 4 mm X 20 cm	
ML-106.43	Ø 4 mm X 25 cm Ø 4 mm X 30 cm	
ML-106.44 ML-106.45	Ø 4 mm X 30 cm	
ML-106.45	Ø 5 mm X 15 cm	
ML-106.47	Ø 5 mm X 20 cm	
ML-106.48	Ø 5 mm X 25 cm	
ML-106.49	Ø 5 mm X 30 cm	
ML-106.50	Ø 5 mm X 35 cm	
ML-106.51	Ø 6 mm X 15 cm	
ML-106.52	Ø 6 mm X 20 cm	
ML-106.53	Ø 6 mm X 25 cm	
ML-106.54	Ø 6 mm X 30 cm	
ML-106.55	Ø 6 mm X 35 cm	



### **INFILTRATION CANNULAS**



Tumescent Infiltrator Angled
This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

	60cc Tommey hub	
ML-101.21	Ø 2.0 mm X 15 cm	
ML-101.22	Ø 2.5 mm X 20 cm	
ML-101.23	Ø 3.0 mm X 20 cm	
ML-101.24	Ø 3.0 mm X 25 cm	

	Treaded Hub	
ML-102.21	Ø 2.0mm X 15 cm	
ML-102.22	Ø 3.0mm X 15 cm	
ML-102.23	Ø 3.0mm X 20 cm	
ML-102.24	Ø 3.0mm X 30 cm	
ML-102.25	Ø 4.0mm X 20 cm	
ML-102.26	Ø 4.0mm X 30 cm	

	Super Luer Lock	
ML-103.21 ML-103.22	Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm	
ML-103.23 ML-103.24	Ø 1.5 mm X 7 cm Ø 1.5 mm X 9 cm	

	Luer lock
	Lasi isok
ML-104.31	Ø 1.0 mm X 4 cm
ML-104.32	Ø 1.0 mm X 5 cm
ML-104.33	Ø 1.0 mm X 7 cm
ML-104.34	Ø 1.5 mm X 7 cm
ML-104.35	Ø 1.5 mm X 9 cm
ML-104.36	Ø 1.5 mm X 10 cm
ML-104.37	Ø 1.5 mm X 12 cm
ML-104.38	Ø 2.0 mm X 10 cm
ML-104.39	Ø 2.0 mm X 12 cm
ML-104.40	Ø 2.0 mm X 15 cm
ML-104.41	Ø 2.0 mm X 20 cm
ML-104.42	Ø 2.5 mm X 15 cm
ML-104.43	Ø 2.5 mm X 20 cm
ML-104.44	Ø 3.0 mm X 15 cm
ML-104.45	Ø 3.0 mm X 20 cm

	Luer lock (Stainless Steel)
ML-105.41	Ø 1.0 mm X 4 cm
ML-105.42	Ø 1.0 mm X 5 cm
ML-105.43	Ø 1.0 mm X 7 cm
ML-105.44	Ø 1.5 mm X 7 cm
ML-105.45	Ø 1.5 mm X 9 cm
ML-105.46	Ø 1.5 mm X 10 cm
ML-105.47	Ø 1.5 mm X 12 cm
ML-105.48	Ø 2.0 mm X 10 cm
ML-105.50	Ø 2.0 mm X 15 cm
ML-105.51	Ø 2.0 mm X 20 cm
ML-105.52	Ø 2.5 mm X 15 cm
ML-105.53	Ø 2.5 mm X 20 cm
ML-105.54	Ø 3.0 mm X 15 cm
ML-105.55	Ø 3.0 mm X 20 cm

	Power Cannula
ML-106.61	Ø 2 mm X 15 cm
ML-106.62	Ø 2 mm X 20 cm
ML-106.63	Ø 2 mm X 25 cm
ML-106.64	Ø 2 mm X 30 cm
ML-106.65	Ø 2 mm X 35 cm
ML-106.66	Ø 3 mm X 15 cm
ML-106.67	Ø 3 mm X 20 cm
ML-106.68	Ø 3 mm X 25 cm
ML-106.69	Ø 3 mm X 30 cm
ML-106.70	Ø 3 mm X 35 cm
ML-106.71	Ø 4 mm X 15 cm
ML-106.72	Ø 4 mm X 20 cm
ML-106.73	Ø 4 mm X 25 cm
ML-106.74	Ø 4 mm X 30 cm
ML-106.75	Ø 4 mm X 35 cm
ML-106.76	Ø 5 mm X 15 cm
ML-106.77	Ø 5 mm X 20 cm
ML-106.78	Ø 5 mm X 25 cm Ø 5 mm X 30 cm
ML-106.79 ML-106.80	Ø 5 mm X 35 cm
ML-106.80	Ø 6 mm X 15 cm
ML-106.82	Ø 6 mm X 20 cm
ML-106.82	Ø 6 mm X 25 cm
ML-106.84	Ø 6 mm X 30 cm
ML-106.85	Ø 6 mm X 35 cm
IVIL-100.00	DO HIIII X 33 GH



### **INFILTRATION CANNULAS**

**Facelift Infiltration Cannula**This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

	60cc Tommey hub	
ML-101.31	Ø 2.0 mm X 15 cm	
ML-101.32	Ø 2.5 mm X 20 cm	
ML-101.33	Ø 3.0 mm X 20 cm	
ML-101.34	Ø 3.0 mm X 25 cm	

	Treaded Hub	
ML-102.31	Ø 2.0mm X 15 cm	
ML-102.32	Ø 3.0mm X 15 cm	
ML-102.33	Ø 3.0mm X 20 cm	
ML-102.34	Ø 3.0mm X 30 cm	
ML-102.35	Ø 4.0mm X 20 cm	
ML-102.36	Ø 4.0mm X 30 cm	

	Super Luer Lock	
ML-103.31 ML-103.32	Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm	
ML-103.33 ML-103.34	Ø 1.5 mm X 7 cm Ø 1.5 mm X 9 cm	

	Luer lock	
ML-104.51	Ø 1.0 mm X 4 cm	
ML-104.52	Ø 1.0 mm X 5 cm	
ML-104.53	Ø 1.0 mm X 7 cm	
ML-104.54	Ø 1.5 mm X 7 cm	
ML-104.55	Ø 1.5 mm X 9 cm	
ML-104.56	Ø 1.5 mm X 10 cm	
ML-104.57	Ø 1.5 mm X 12 cm	
ML-104.58	Ø 2.0 mm X 10 cm	
ML-104.59	Ø 2.0 mm X 12 cm	
ML-104.50	Ø 2.0 mm X 15 cm	
ML-104.51	Ø 2.0 mm X 20 cm	
ML-104.52	Ø 2.5 mm X 15 cm	
ML-104.53	Ø 2.5 mm X 20 cm	
ML-104.54	Ø 3.0 mm X 15 cm	
ML-104.55	Ø 3.0 mm X 20 cm	

	Luer lock (Stainless Steel)
ML-105.61	Ø 1.0 mm X 4 cm
ML-105.62	Ø 1.0 mm X 5 cm
ML-105.63	Ø 1.0 mm X 7 cm
ML-105.64	Ø 1.5 mm X 7 cm
ML-105.65	Ø 1.5 mm X 9 cm
ML-105.66	Ø 1.5 mm X 10 cm
ML-105.67	Ø 1.5 mm X 12 cm
ML-105.68	Ø 2.0 mm X 10 cm
ML-105.69	Ø 2.0 mm X 12 cm
ML-105.70	Ø 2.0 mm X 15 cm
ML-105.71	Ø 2.0 mm X 20 cm
ML-105.72	Ø 2.5 mm X 15 cm
ML-105.73	Ø 2.5 mm X 20 cm
ML-105.74	Ø 3.0 mm X 15 cm
ML-105.75	Ø 3.0 mm X 20 cm

	Power Cannula
ML-106.91	Ø 2 mm X 15 cm
ML-106.92	Ø 2 mm X 20 cm
ML-106.93	Ø 2 mm X 25 cm
ML-106.94	Ø 2 mm X 30 cm
ML-106.95	Ø 2 mm X 35 cm
ML-106.96	Ø 3 mm X 15 cm
ML-106.97	Ø 3 mm X 20 cm
ML-106.98	Ø 3 mm X 25 cm
ML-106.99	Ø 3 mm X 30 cm
ML-106.100	Ø 3 mm X 35 cm
ML-106.101	Ø 4 mm X 15 cm
ML-106.102	Ø 4 mm X 20 cm
ML-106.103	Ø 4 mm X 25 cm
ML-106.104	Ø 4 mm X 30 cm
ML-106.105 ML-106.106	Ø 4 mm X 35 cm Ø 5 mm X 15 cm
ML-106.107	Ø 5 mm X 20 cm
ML-106.108	Ø 5 mm X 25 cm
ML-106.109	Ø 5 mm X 30 cm
ML-106.110	Ø 5 mm X 35 cm
ML-106.111	Ø 6 mm X 15 cm
ML-106.112	Ø 6 mm X 20 cm
ML-106.113	Ø 6 mm X 25 cm
ML-106.114	Ø 6 mm X 30 cm
ML-106.115	Ø 6 mm X 35 cm


medastepro@gmail.com





General suction cannula are ideal tool for large volume Harvesting with ease perfection. These cannulas are designed to facilitate you to get perfectly designed goals in minimum time.





### Stevens Style Cannula

	60cc Tommey hub	
ML-201.01 ML-201.02 ML-201.03 ML-201.04	Ø 2.0 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	

	Treaded Hub	
ML-202.01	Ø 3.0 mm X 20 cm	
ML-202.02	Ø 3.0 mm X 30 cm	
ML-202.03	Ø 4.0 mm X 30 cm	
ML-202.04	Ø 5.0 mm X 30 cm	
ML-202.05	Ø 6.0 mm X 30 cm	

	Super Luer Lock
ML-203.01	Ø 1.0 mm X 5 cm
ML-203.02	Ø 1.0 mm X 7 cm
ML-203.03	Ø 1.5 mm X 7 cm
ML-203.04	Ø 1.5 mm X 9 cm
ML-203.05	Ø 1.5 mm X 7 cm

	Luer lock
ML-204.01 ML-204.02 ML-204.03 ML-204.04 ML-204.05	Ø 1.0 mm X 4 cm Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm Ø 1.5 mm X 7 cm Ø 1.5 mm X 9 cm
ML-204.06 ML-204.07 ML-204.08	Ø 1.5 mm X 10 cm Ø 1.5 mm X 12 cm Ø 2.0 mm X 10 cm
ML-204.09 ML-204.10 ML-204.11 ML-204.12	Ø 2.0 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.0 mm X 20 cm Ø 2.5 mm X 15 cm
ML-204.13 ML-204.14 ML-204.15	Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)		
ML-205.01 ML-205.02 ML-205.03 ML-205.04 ML-205.05 ML-205.06 MI-205.07	Ø 1.0 mm X 4 cm Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm Ø 1.5 mm X 7 cm Ø 1.5 mm X 9 cm Ø 1.5 mm X 10 cm Ø 1.5 mm X 12 cm	
ML-205.08 ML-205.09 ML-205.10 ML-205.11 ML-205.12 ML-205.13 ML-205.14 ML-205.15	Ø 2.0 mm X 10 cm Ø 2.0 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.0 mm X 20 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm	

	Power Cannula
ML-206.01	Ø 2 mm X 15 cm
ML-206.02	Ø 2 mm X 20 cm
ML-206.03	Ø 2 mm X 25 cm
ML-206.04	Ø 2 mm X 30 cm
ML-206.05	Ø 2 mm X 35 cm
ML-206.06	Ø 3 mm X 15 cm
ML-206.07	Ø 3 mm X 20 cm
ML-206.08	Ø 3 mm X 25 cm
ML-206.09	Ø 3 mm X 30 cm
ML-206.10	Ø 3 mm X 35 cm
ML-206.11	Ø 4 mm X 15 cm
ML-206.12	Ø 4 mm X 20 cm
ML-206.13	Ø 4 mm X 25 cm
ML-206.14	Ø 4 mm X 30 cm
ML-206.15	Ø 4 mm X 35 cm
ML-206.16	Ø 5 mm X 15 cm
ML-206.17	Ø 5 mm X 20 cm
ML-206.18	Ø 5 mm X 25 cm
ML-206.19	Ø 5 mm X 30 cm
ML-206.20	Ø 5 mm X 35 cm
ML-206.21	Ø 6 mm X 15 cm
ML-206.22	Ø 6 mm X 20 cm
ML-206.23	Ø 6 mm X 25 cm
ML-206.24	Ø 6 mm X 30 cm
ML-206.25	Ø 6 mm X 35 cm





### Cobra Bibevel Style

More aggressive than the Standard, Pyramid or Mercedes, this cannula has a tapered tip with two ports and works effectively within dense fibrous tissue.

	60cc Tommey hub		
ML-201.11	Ø 2.5 mm X 15 cm5		
ML-201.12	Ø 2.5 mm X 20 cm		
ML-201.13	Ø 3.0 mm X 15 cm		
ML-201.14	Ø 3.0 mm X 20 cm		
ML-201.15	Ø 3.0 mm X 25 cm		
ML-201.16	Ø 4.0 mm X 25 cm		
ML-201.17	Ø 4.0 mm X 30 cm		
ML-201.18	Ø 4.0 mm X 35 cm		
ML-201.19	Ø 4.5 mm X 30 cm		
ML-201.20	Ø 4.5 mm X 35 cm		

Treaded Hub		
ML-202.11 ML-202.12 ML-202.13	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm	
ML-202.14	Ø 4.0 mm X 25 cm	

MI 000 44	Super Luer Lock		
ML-203.11 Ø 2.0 mm X 15 cm ML-203.12 Ø 2.5 mm X 15 cm ML-203.13 Ø 2.5 mm X 20 cm ML-203.14 Ø 3.0 mm X 15 cm ML-203.15 Ø 3.0 mm X 20 cm ML-203.16 Ø 3.0 mm X 25 cm	mm X 20 cm mm X 15 cm mm X 20 cm	ML-203.13 ML-203.14 ML-203.15	

Luer lock		
ML-204.21	Ø 2.0 mm X 15 cm	
ML-204.22	Ø 2.5 mm X 15 cm	
ML-204.23	Ø 2.5 mm X 20 cm	
ML-204.24	Ø 3.0 mm X 15 cm	
ML-204.25	Ø 3.0 mm X 20 cm	
ML-204.26	Ø 3.0 mm X 25 cm	

Luer lock	(Stainless Steel)	
ML-205.21 ML-205.22 ML-205.23 ML-205.24 ML-205.25 ML-205.26	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	

Power	Cannula	
ML-206.31	1 Ø 2mm X 15 cm	
ML-206.32	2 Ø 2mm X 20 cm	
ML-206.33	8 Ø 2mm X 25 cm	
ML-206.34	Ø 2mm X 30 cm	
ML-206.35	Ø 2mm X 35 cm	
ML-206.36	Ø 3mm X 15 cm	
ML-206.37	7 Ø 3mm X 20 cm	
ML-206.38	3 Ø 3mm X 25 cm	
ML-206.39	9 Ø 3mm X 30 cm	
ML-206.40		
ML-206.41		
ML-206.42		
ML-206.43		
ML-206.44		
ML-206.45		
ML-206.46		
ML-206.47		
ML-206.48		
ML-206.49		
ML-206.50		
ML-206.51		
ML-206.52		
ML-206.53		
ML-206.54		
ML-206.55	Ø 6mm X 35 cm	





**Cobra Style Round Tip**Less aggressive than the Cobra Bibevel this cell friendly cannula has a rounded tip with two ports and works effectively within dense fibrous tissue during liposuction procedures.

60cc Tommey hub		
ML-201.21 ML-201.22 ML-201.23 ML-201.24	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-201.25 ML-201.26 ML-201.27 ML-201.28 ML-201.29 ML-201.30	Ø 3.0 mm X 25 cm Ø 4.0 mm X 25 cm Ø 4.0 mm X 30 cm Ø 4.0 mm X 35 cm Ø 4.5 mm X 30 cm Ø 4.5 mm X 35 cm	

Treaded Hub		
ML-202.21 ML-202.22 ML-202.23 ML-202.24	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	


Super Luer Lock		
ML-203.21	Ø 2.0 mm X 15 cm	
ML-203.22	Ø 2.5 mm X 15 cm	
ML-203.23	Ø 2.5 mm X 20 cm	
ML-203.24	Ø 3.0 mm X 15 cm	
ML-203.25	Ø 3.0 mm X 20 cm	
ML-203.26	Ø 3.0 mm X 25 cm	

	Luer lock	
ML-204.31	Ø 2.0 mm X 15 cm	
ML-204.32 ML-204.33	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm	
ML-204.34	Ø 3.0 mm X 15 cm	
ML-204.35 ML-204.36	Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	
WIL-204.30	Ø 3.0 IIIII A 23 GII	

Luer lock (Stainless Steel)		
ML-205.31 ML-205.32 ML-205.33 ML-205.34 ML-205.35	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-205.36	Ø 3.0 mm X 25 cm	

Power Ca	nnula	
ML-206.61	Ø 2 mm X 15 cm	
ML-206.62	Ø 2 mm X 20 cm	
ML-206.63	Ø 2 mm X 25 cm	
ML-206.64	Ø 2 mm X 30 cm	
ML-206.65	Ø 2 mm X 35 cm	
ML-206.66	Ø 3 mm X 15 cm	
ML-206.67	Ø 3 mm X 20 cm	
ML-206.68	Ø 3 mm X 25 cm	
ML-206.69	Ø 3 mm X 30 cm	
ML-206.70	Ø 3 mm X 35 cm	
ML-206.71	Ø 4 mm X 15 cm	
ML-206.72	Ø 4 mm X 20 cm	
ML-206.73	Ø 4 mm X 25 cm	
ML-206.74	Ø 4 mm X 30 cm	
ML-206.75	Ø 4 mm X 35 cm	
ML-206.76	Ø 5 mm X 15 cm	
ML-206.77	Ø 5 mm X 20 cm	
ML-206.78	Ø 5 mm X 25 cm	
ML-206.79	Ø 5 mm X 30 cm	
ML-206.80	Ø 5 mm X 35 cm	
ML-206.81	Ø 6 mm X 15 cm	
ML-206.82	Ø 6 mm X 20 cm	
ML-206.83	Ø 6 mm X 25 cm	
ML-206.84	Ø 6 mm X 30 cm	
ML-206.85	Ø 6 mm X 35 cm	







**Standard, Single Hole**This is the most common cell friendly cannula design for soft tissue removal. The single port is located 180 degrees from the thumb rest on the hub.

60cc Tommey hub		
ML-201.31	Ø 2.5 mm X 15 cm	
ML-201.32	Ø 2.5 mm X 20 cm	
ML-201.33	Ø 3.0 mm X 15 cm	
ML-201.34	Ø 3.0 mm X 20 cm	
ML-201.35	Ø 3.0 mm X 25 cm	
ML-201.36	Ø 4.0 mm X 25 cm	
ML-201.37	Ø 4.0 mm X 30 cm	
ML-201.38	Ø 4.0 mm X 35 cm	
ML-201.39	Ø 4.5 mm X 30 cm	
ML-201.40	Ø 4.5 mm X 35 cm	

Treaded Hub		
ML-202.31 ML-202.32 ML-202.33 ML-202.34	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

Super Luer Lock		
ML-203.32	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.0 mm X 20 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-203.38	Ø 3.0 mm X 25 cm	

	Luer lock	
ML-204.41	Ø 1.5 mm X 12 cm	
ML-204.42	Ø 2.0 mm X 15 cm	
ML-204.43	Ø 2.0 mm X 20 cm	
ML-204.44	Ø 2.5 mm X 15 cm	
ML-204.45	Ø 2.5 mm X 20 cm	
ML-204.46	Ø 3.0 mm X 15 cm	
ML-204.47	Ø 3.0 mm X 20 cm	
ML <b>-</b> 204.48	Ø 3.0 mm X 25 cm	

Luer lock (Stainless Steel)		
ML-205.41 ML-205.42 ML-205.43 ML-205.44 ML-205.45 ML-205.46 ML-205.47	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.0 mm X 20 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-205.48	Ø 3.0 mm X 25 cm	

	Power Cannula
ML-206.91	Ø 2 mm X 15 cm
ML-206.92	Ø 2 mm X 20 cm
ML-206.93	Ø 2 mm X 25 cm
ML-206.94	Ø 2 mm X 30 cm
ML-206.95	Ø 2 mm X 35 cm
ML-206.96	Ø 3 mm X 15 cm
ML-206.97	Ø 3 mm X 20 cm
ML-206.98	Ø 3 mm X 25 cm
ML-206.99	Ø 3 mm X 30 cm
ML-206.100	Ø 3 mm X 35 cm
ML-206.101	Ø 4 mm X 15 cm
ML-206.102	Ø 4 mm X 20 cm
ML-206.103	Ø 4 mm X 25 cm
ML-206.104	Ø 4 mm X 30 cm
ML-206.105	Ø 4 mm X 35 cm
ML-206.106	Ø 5 mm X 15 cm
ML-206.107	Ø 5 mm X 20 cm
ML-206.108	Ø 5 mm X 25 cm
ML-206.109	Ø 5 mm X 30 cm
ML-206.110	Ø 5 mm X 35 cm
ML-206.111	Ø 6 mm X 15 cm
ML-206.112	Ø 6 mm X 20 cm
ML-206.113	Ø 6 mm X 25 cm
ML-206.114	Ø 6 mm X 30 cm
ML-206.115	Ø 6 mm X 35 cm







**Pyramid Holes Cannula**This commonly used design has three ports located 180 degrees from the thumb print. It is a slightly more aggressive cannula than the Standard during liposuction procedures.

60cc Tommey hub		
ML-201.41	Ø 2.5 mm X 15 cm	
ML-201.42	Ø 2.5 mm X 20 cm	
ML-201.43	Ø 3.0 mm X 15 cm	
ML-201.44	Ø 3.0 mm X 20 cm	
ML-201.45	Ø 3.0 mm X 25 cm	
ML-201.46	Ø 4.0 mm X 25 cm	
ML-201.47	Ø 4.0 mm X 30 cm	
ML-201.48	Ø 4.0 mm X 35 cm	
ML-201.49	Ø 4.5 mm X 30 cm	
ML-201.50	Ø 4.5 mm X 35 cm	

Treaded Hub		
ML-202.41 ML-202.42 ML-202.43 ML-202.44	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	


Super Luer Lock		
ML-203.41 ML-203.42 ML-203.43 ML-203.44 ML-203.45 ML-203.46	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-203.47	Ø 3.0 mm X 25 cm	

	Luer lock	
ML-204.51 ML-204.52 ML-204.53	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm	
ML-204.53 ML-204.54 ML-204.55	Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm	
ML-204.56 ML-204.57	Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	

	Luer lock (Stainless Steel)
ML-205.51	Ø 1.5 mm X 12 cm
ML-205.52	Ø 2.0 mm X 15 cm
ML-205.53	Ø 2.0 mm X 20 cm
ML-205.54	Ø 2.5 mm X 15 cm
ML-205.55	Ø 2.5 mm X 20 cm
ML-205.56	Ø 3.0 mm X 15 cm
ML-205.57	Ø 3.0 mm X 20 cm
ML-205.58	Ø 3.0 mm X 25 cm

	Power Cannula
ML-206.121	Ø 2 mm X 15 cm
ML-206.122	Ø 2 mm X 20 cm
ML-206.123	Ø 2 mm X 25 cm
ML-206.124	Ø 2 mm X 30 cm
ML-206.125	Ø 2 mm X 35 cm
ML-206.126	Ø 3 mm X 15 cm
ML-206.127	Ø 3 mm X 20 cm
ML-206.128	Ø 3 mm X 25 cm
ML-206.129	Ø 3 mm X 30 cm
ML-206.130	Ø 3 mm X 35 cm
ML-206.131	Ø 4 mm X 15 cm
ML-206.132	Ø 4 mm X 20 cm
ML-206.133	Ø 4 mm X 25 cm
ML-206.134	Ø 4 mm X 30 cm
ML-206.135	Ø 4 mm X 35 cm
ML-206.136	Ø 5 mm X 15 cm
ML-206.137	Ø 5 mm X 20 cm
ML-206.138	Ø 5 mm X 25 cm
ML-206.139	Ø 5 mm X 30 cm
ML-206.140	Ø 5 mm X 35 cm
ML-206.141	Ø 6 mm X 15 cm
ML-206.142	Ø 6 mm X 20 cm
ML-206.143	Ø 6 mm X 25 cm
ML-206.144	Ø 6 mm X 30 cm
ML-206.145	Ø 6 mm X 35 cm





#### Mercedes

This circumferential three port configuration cell friendly liposuction cannula is effective for multi-plane soft tissue harvesting.

	60cc Tommey hub	
ML-201.51	Ø 2.5 mm X 15 cm	
ML-201.52	Ø 2.5 mm X 20 cm	
ML-201.53	Ø 3.0 mm X 15 cm	
ML-201.54	Ø 3.0 mm X 20 cm	
ML-201.55	Ø 3.0 mm X 25 cm	
ML-201.56	Ø 4.0 mm X 25 cm	
ML-201.57	Ø 4.0 mm X 30 cm	
ML-201.58	Ø 4.0 mm X 35 cm	
ML-201.59	Ø 4.5 mm X 30 cm	
ML-201.60	Ø 4.5 mm X 35 cm	

	Treaded Hub	
ML-202.51 ML-202.52 ML-202.53 ML-202.54	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

	Super Luer Lock	
ML-203.51 ML-203.52 ML-203.53 ML-203.54 ML-203.55 ML-203.56 ML-203.57	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	

	Luer lock	
ML-204.61 ML-204.62 ML-204.63 ML-204.64 ML-204.65 ML-204.66	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-204.67	Ø 3.0 mm X 25 cm	

	Luer lock (Stainless Steel)
ML-205.61	Ø 1.5 mm X 12 cm
ML-205.62	Ø 2.0 mm X 15 cm
ML-205.63	Ø 2.0 mm X 20 cm
ML-205.64	Ø 2.5 mm X 15 cm
ML-205.65	Ø 2.5 mm X 20 cm
ML-205.66	Ø 3.0 mm X 15 cm
ML-205.67	Ø 3.0 mm X 20 cm
ML-205.68	Ø 3.0 mm X 25 cm

	Power Cannula
ML-206.151	Ø 2 mm X 15 cm
ML-206.152	Ø 2 mm X 20 cm
ML-206.153	Ø 2 mm X 25 cm
ML-206.154	Ø 2 mm X 30 cm
ML-206.155	Ø 2 mm X 35 cm
ML-206.156	Ø 3 mm X 15 cm
ML-206.157	Ø 3 mm X 20 cm
ML-206.158	Ø 3 mm X 25 cm
ML-206.159	Ø 3 mm X 30 cm
ML-206.160	Ø 3 mm X 35 cm
ML-206.161	Ø 4 mm X 15 cm
ML-206.162	Ø 4 mm X 20 cm
ML-206.163	Ø 4 mm X 25 cm
ML-206.164	Ø 4 mm X 30 cm
ML-206.165	Ø 4 mm X 35 cm
ML-206.166	Ø 5 mm X 15 cm
ML-206.167	Ø 5 mm X 20 cm
ML-206.168	Ø 5 mm X 25 cm
ML-206.169	Ø 5 mm X 30 cm
ML-206.170	Ø 5 mm X 35 cm
ML-206.171	Ø 6 mm X 15 cm
ML-206.172	Ø 6 mm X 20 cm
ML-206.173	Ø 6 mm X 25 cm
ML-206.174	Ø 6 mm X 30 cm
ML-206.175	Ø 6 mm X 35 cm





3 Spiral Port Cannula (Carraway Style)
Spiral three port design on this Cell Friendly liposuction cannula is effective for multi-plane adipose tissue harvesting.

	60cc Tommey hub	
ML-201.61 ML-201.62 ML-201.63 ML-201.65 ML-201.65 ML-201.66 ML-201.68 ML-201.68	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 25 cm Ø 4.0 mm X 30 cm Ø 4.0 mm X 30 cm	
ML-201.70	Ø 4.5 mm X 35 cm  Treaded Hub	
ML-202.61 ML-202.62 ML-202.63 ML-202.64	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

	Super Luer Lock	
ML-203.61 ML-203.62 ML-203.63 ML-203.64 ML-203.65	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm	
ML-203.66 ML-203.67	Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	


	Luer lock	
ML-204.71 ML-204.72 ML-204.73 ML-204.74	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm	
ML-204.75 ML-204.76 ML-204.77	Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	

	Luer lock (Stainless Steel)
ML-205.71	Ø 1.5 mm X 12 cm
ML-205.72	Ø 2.0 mm X 15 cm
ML-205.73	Ø 2.0 mm X 20 cm
ML-205.74	Ø 2.5 mm X 15 cm
ML-205.75	Ø 2.5 mm X 20 cm
ML-205.76	Ø 3.0 mm X 15 cm
ML-205.77	Ø 3.0 mm X 20 cm
ML-205.77 ML-205.78	Ø 3.0 mm X 25 cm

Power Cannula			
ML-206.181	Ø 2 mm X 15 cm		
ML-206.182	Ø 2 mm X 20 cm		
ML-206.183	Ø 2 mm X 25 cm		
ML-206.184	Ø 2 mm X 30 cm		
ML-206.185	Ø 2 mm X 35 cm		
ML-206.186	Ø 3 mm X 15 cm		
ML-206.187	Ø 3 mm X 20 cm		
ML-206.188	Ø 3 mm X 25 cm		
ML-206.189	Ø 3 mm X 30 cm		
ML-206.190	Ø 3 mm X 35 cm		
ML-206.191	Ø 4 mm X 15 cm		
ML-206.192	Ø 4 mm X 20 cm		
ML-206.192 ML-206.193 ML-206.194 ML-206.195 ML-206.196 ML-206.197	Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm		
ML-206.198	Ø 5 mm X 25 cm		
ML-206.199	Ø 5 mm X 30 cm		
ML-206.200	Ø 5 mm X 35 cm		
ML-206.201	Ø 6 mm X 15 cm		
ML-206.202	Ø 6 mm X 20 cm		
ML-206.203	Ø 6 mm X 25 cm		
ML-206.204	Ø 6 mm X 30 cm		
ML-206.205	Ø 6 mm X 35 cm		



### MEDASTE VOTRE APPRENTISSACE NOTRE DEVOIR GENERAL SUCTION CANNULAS





Six Spiral Port (Adjust below 3 Spiral Ports)
Six port spiral configuration on this Cell Friendly liposuction cannula harvests maximum tissue volume in minimal time.

6	60cc smart hub		
ML	-201.91	Ø 2.5 mm X 15 cm	
ML	-201.92	Ø 2.5 mm X 20 cm	
ML	-201.93	Ø 3.0 mm X 15 cm	
ML	-201.94	Ø 3.0 mm X 25 cm	
ML	-201.95	Ø 4.0 mm X 25 cm	

Treaded Hub			
ML-202.81 ML-202.82 ML-202.83 ML-202.84	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm		

	Super Luer Lock	
ML-203.81	Ø 2.5 mm X 15 cm	
ML-203.82	Ø 3.0 mm X 15 cm	
ML-203.83	Ø 3.0 mm X 25 cm	
ML-203.84	Ø 4.0 mm X 15 cm	
ML-203.85	Ø 4.0 mm X 25 cm	


Luer lock	
ML-204.91	Ø 2.5 mm X 15 cm
ML-204.92	Ø 3.0 mm X 15 cm
ML-204.93	Ø 3.0 mm X 25 cm
ML-204.94	Ø 4.0 mm X 15 cm
ML-204.95	Ø 4.0 mm X 25 cm

Luer lock (Stainless Steel)			
ML-205.231	Ø 2.5 mm X 15 cm		
ML-205.232 ML-205.233	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm		
ML-205.234	Ø 4.0 mm X 15 cm		
ML-205.235	Ø 4.0 mm X 25 cm		

Power Ca	nnula
ML-206.901	2mm OD X 15 cm
ML-206.902	2mm OD X 20 cm
ML-206.903	2mm OD X 25 cm
ML-206.904	2mm OD X 30 cm
ML-206.905	2mm OD X 35 cm
ML-206.906	3mm OD X 15 cm
ML-206.907	3mm OD X 20 cm
ML-206.908	3mm OD X 25 cm
ML-206.909	3mm OD X 30 cm
ML-206.910	3mm OD X 35 cm
ML-206.911	4mm OD X 15 cm
ML-206.912	4mm OD X 20 cm
ML-206.913	4mm OD X 25 cm
ML-206.914	4mm OD X 30 cm
ML-206.915	4mm OD X 35 cm
ML-206.916	5mm OD X 15 cm
ML-206.917	5mm OD X 20 cm
ML-206.918	5mm OD X 25 cm
ML-206.919	5mm OD X 30 cm
ML-206.920	5mm OD X 35 cm
ML-206.921	6mm OD X 15 cm
ML-206.922	6mm OD X 20 cm
ML-206.923	6mm OD X 25 cm
ML-206.924	6mm OD X 30 cm
ML-206.925	6mm OD X 35 cm

# MEDASTE VOTRE APPRENTISSAGE, NOTRE DEVOIR

## GENERAL SUCTION CANNULAS





#### Three Port Cannula

This body sculpting instrument is excellent for harvesting small aliquots of tissue during liposuction procedures.

	60cc Tommey hub	
ML-201.71	Ø 2.5 mm X 15 cm	
ML-201.72	Ø 2.5 mm X 20 cm	
ML-201.73	Ø 3.0 mm X 15 cm	
ML-201.74	Ø 3.0 mm X 20 cm	
ML-201.75	Ø 3.0 mm X 25 cm	
ML-201.76	Ø 4.0 mm X 25 cm	
ML-201.77	Ø 4.0 mm X 30 cm	
ML-201.78	Ø 4.0 mm X 35 cm	
ML-201.79	Ø 4.5 mm X 30 cm	
ML-201.80	Ø 4.5 mm X 35 cm	

Luer lock (Stainless Steel)			
ML-205.81	Ø 1.5 mm X 12 cm		
ML-205.82	Ø 2.0 mm X 15 cm		
ML-205.83	Ø 2.0 mm X 20 cm		
ML-205.84	Ø 2.5 mm X 15 cm		
ML-205.85	Ø 2.5 mm X 20 cm		
ML-205.86	Ø 3.0 mm X 15 cm		
ML-205.87	Ø 3.0 mm X 20 cm		
ML-205.88	Ø 3.0 mm X 25 cm		

Treaded Hub		
ML-202.71 ML-202.72 ML-202.73 ML-202.74	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

	Super Luer Lock	
ML-203.71 ML-203.72 ML-203.73 ML-203.74 ML-203.75 ML-203.76	Ø 1.5 mm X 12 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-203.77	Ø 3.0 mm X 25 cm	

	Power Cannula
ML-206.211	Ø 2 mm X 15 cm
ML-206.212	Ø 2 mm X 20 cm
ML-206.213	Ø 2 mm X 25 cm
ML-206.214	Ø 2 mm X 30 cm
ML-206.215	Ø 2 mm X 35 cm
ML-206.216	Ø 3 mm X 15 cm
ML-206.217	Ø 3 mm X 20 cm
ML-206.218	Ø 3 mm X 25 cm
ML-206.219	Ø 3 mm X 30 cm
ML-206.220	Ø 3 mm X 35 cm
ML-206.221	Ø 4 mm X 15 cm
ML-206.222	Ø 4 mm X 20 cm
ML-206.223	Ø 4 mm X 25 cm
ML-206.224	Ø 4 mm X 30 cm
ML-206.225	Ø 4 mm X 35 cm
ML-206.226	Ø 5 mm X 15 cm
ML-206.227	Ø 5 mm X 20 cm
ML-206.228	Ø 5 mm X 25 cm
ML-206.229	Ø 5 mm X 30 cm
ML-206.230	Ø 5 mm X 35 cm
ML-206.231	Ø 6 mm X 15 cm
ML-206.232	Ø 6 mm X 20 cm
ML-206.233	Ø 6 mm X 25 cm
ML-206.234	Ø 6 mm X 30 cm
ML-206.235	Ø 6 mm X 35 cm





#### **Basket Shaped Cannula**

The Cell Friendly Basket Cannula is a versatile liposculpting tool. It can operate as an equalizer without vacuum to treat adipose tissue and skin irregularities in revision liposuction and fat transfer procedures. The Basket Cannula is also adept under vacuum as a fat harvesting cannula. For finishing, this cannula can be used as an injector for equalizing and adding adipose tissue to sculpt smoother contures.

60cc Tommey hub					
ML-201.81	Ø 3.0 mm X 15 cm				
ML-201.82	Ø 3.0 mm X 20 cm				
ML-201.83	Ø 3.0 mm X 25 cm				
ML-201.84	Ø 4.0 mm X 25 cm				
ML-201.85	Ø 4.0 mm X 30 cm				
ML-201.86	Ø 4.0 mm X 35 cm				



### Hole on Tip-Style I


Power Cannula							
ML-206.251 ML-206.252 ML-206.253 ML-206.254 ML-206.255 ML-206.256 ML-206.257 ML-206.258 ML-206.259	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.260 ML-206.261 ML-206.262 ML-206.263 ML-206.264 ML-206.265 ML-206.266 ML-206.267	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.268 ML-206.269 ML-206.270 ML-206.271 ML-206.272 ML-206.273 ML-206.274 ML-206.275	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		



#### Hole on tip-Style II

Power Cannula							
ML-206.281 ML-206.282 ML-206.283 ML-206.284 ML-206.285 ML-206.286 ML-206.287 ML-206.288 ML-206.289	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.290 ML-206.291 ML-206.292 ML-206.293 ML-206.294 ML-206.295 ML-206.296 ML-206.297	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.298 ML-206.299 ML-206.300 ML-206.301 ML-206.302 ML-206.303 ML-206.304 ML-206.305	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		







### **Deep Sharp Cut Dissector**

Power Cannula							
ML-206.311 ML-206.312 ML-206.313 ML-206.314 ML-206.315 ML-206.316 ML-206.317 ML-206.318 ML-206.319	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.320 ML-206.321 ML-206.322 ML-206.323 ML-206.324 ML-206.325 ML-206.326 ML-206.327	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.328 ML-206.329 ML-206.330 ML-206.331 ML-206.332 ML-206.333 ML-206.334 ML-206.335	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		




### **Bi-Port Cannula**

Luer lock (Stainless Steel)								
	10cm	15cm	20cm	25cm				
2.0 Ø 2.5 Ø	ML-205.91 ML-205.92	ML-205.101 ML-205.102	ML-205.111 ML-205.112					
3.0 Ø	ML-205.93	ML-205.103	ML-205.113					
3.5 Ø	ML-205.94	ML-205.104	ML-205.114	ML-205.124				
4.0 Ø	ML-205.95	ML-205.105	ML-205.115	ML-205.125				

Power Cannula							
ML-206.351 ML-206.352 ML-206.353 ML-206.354 ML-206.355 ML-206.356 ML-206.357 ML-206.358 ML-206.359	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.360 ML-206.361 ML-206.362 ML-206.363 ML-206.364 ML-206.365 ML-206.366 ML-206.367	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.368 ML-206.369 ML-206.370 ML-206.371 ML-206.372 ML-206.373 ML-206.374 ML-206.375	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		





### Tappered Tip Single Central Hole

Power Cannula							
ML-206.391 ML-206.392 ML-206.393 ML-206.394 ML-206.395 ML-206.396 ML-206.397 ML-206.398 ML-206.399	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.400 ML-206.401 ML-206.402 ML-206.403 ML-206.404 ML-206.405 ML-206.406 ML-206.407	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.408 ML-206.409 ML-206.410 ML-206.411 ML-206.412 ML-206.413 ML-206.414	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		



### **Tapperd Tip Two Spiral Holes**

Power Cannula							
ML-206.431 ML-206.432 ML-206.433 ML-206.434 ML-206.435 ML-206.436 ML-206.437 ML-206.438 ML-206.439	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.440 ML-206.441 ML-206.442 ML-206.443 ML-206.444 ML-206.445 ML-206.446 ML-206.447	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.448 ML-206.449 ML-206.450 ML-206.451 ML-206.452 ML-206.453 ML-206.454 ML-206.455	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		



Power Cannula							
ML-206.471 ML-206.472 ML-206.473 ML-206.474 ML-206.475 ML-206.476 ML-206.477 ML-206.478 ML-206.479	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.480 ML-206.481 ML-206.482 ML-206.483 ML-206.484 ML-206.485 ML-206.486 ML-206.487	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.488 ML-206.489 ML-206.490 ML-206.491 ML-206.492 ML-206.493 ML-206.494 ML-206.495	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		









Fig. 2Ech


Power Cannula						
ML-206.501 ML-206.502 ML-206.503 ML-206.504 ML-206.505 ML-206.506 ML-206.507 ML-206.508 ML-206.509	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 25 cm	ML-206.510 ML-206.511 ML-206.512 ML-206.513 ML-206.514 ML-206.515 ML-206.516 ML-206.517	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.518 ML-206.519 ML-206.520 ML-206.521 ML-206.522 ML-206.523 ML-206.524 ML-206.525	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm	





Power Cannula							
ML-206.551 ML-206.552 ML-206.553 ML-206.554 ML-206.555 ML-206.556 ML-206.557 ML-206.558 ML-206.559	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.560 ML-206.561 ML-206.562 ML-206.563 ML-206.564 ML-206.565 ML-206.566	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.568 ML-206.569 ML-206.570 ML-206.571 ML-206.572 ML-206.573 ML-206.574 ML-206.575	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm		







#### **One Central Port-Two Parallel Ports**

Power Cannula						
ML-206.591 ML-206.592 ML-206.593 ML-206.594 ML-206.595 ML-206.596 ML-206.597 ML-206.598 ML-206.599	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.600 ML-206.601 ML-206.602 ML-206.603 ML-206.604 ML-206.605 ML-206.606 ML-206.607	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.608 ML-206.609 ML-206.610 ML-206.611 ML-206.612 ML-206.613 ML-206.614 ML-206.615	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm	



### Two Parallel Port- One Central Port

Power Cannula					
ML-206.651 ML-206.652 ML-206.653 ML-206.654 ML-206.655 ML-206.656 ML-206.657 ML-206.658 ML-206.659	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.660 ML-206.661 ML-206.662 ML-206.663 ML-206.664 ML-206.665 ML-206.666	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.668 ML-206.669 ML-206.670 ML-206.671 ML-206.672 ML-206.673 ML-206.674 ML-206.675	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm















Fig. 3Me

Fig. 3Hr

Power Cannula						
ML-206.721 ML-206.722 ML-206.723 ML-206.724 ML-206.725 ML-206.726 ML-206.727 ML-206.728 ML-206.729	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 20 cm Ø 3 mm X 30 cm	ML-206.730 ML-206.731 ML-206.732 ML-206.733 ML-206.734 ML-206.735 ML-206.736 ML-206.737	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.738 ML-206.739 ML-206.740 ML-206.741 ML-206.742 ML-206.743 ML-206.744	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm	

Fig. Riz

	Power Cannula				
ML-206.751	Ø 2 mm X 15 cm	ML-206.760	Ø 3 mm X 35 cm	ML-206.768	Ø 5 mm X 25 cm
ML-206.752	Ø 2 mm X 20 cm	ML-206.761	Ø 4 mm X 15 cm	ML-206.769	Ø 5 mm X 30 cm
ML-206.753	Ø 2 mm X 25 cm	ML-206.762	Ø 4 mm X 20 cm	ML-206.770	Ø 5 mm X 35 cm
ML-206.754	Ø 2 mm X 30 cm	ML-206.763	Ø 4 mm X 25 cm	ML-206.771	Ø 6 mm X 15 cm
ML-206.755	Ø 2 mm X 35 cm	ML-206.764	Ø 4 mm X 30 cm	ML-206.772	Ø 6 mm X 20 cm
ML-206.756	Ø 3 mm X 15 cm	ML-206.765	Ø 4 mm X 35 cm	ML-206.773	Ø 6 mm X 25 cm
ML-206./5/	Ø 3 mm X 20 cm	ML-206.766	Ø 5 mm X 15 cm	ML-206.774	Ø 6 mm X 30 cm
ML-206.758 ML-206.759	Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.767	Ø 5 mm X 20 cm	ML-206.775	Ø 6 mm X 35 cm

medastepro@gmail.com







Fig. 5Ai

Power Cannula						
ML-206.721 ML-206.722 ML-206.723 ML-206.724 ML-206.725 ML-206.726 ML-206.727 ML-206.727	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.730 ML-206.731 ML-206.732 ML-206.733 ML-206.734 ML-206.735 ML-206.736 ML-206.737	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.738 ML-206.739 ML-206.740 ML-206.741 ML-206.742 ML-206.743 ML-206.744 ML-206.745	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm	



Fig. Mr

Power Cannula						
ML-206.751 ML-206.752 ML-206.753 ML-206.754 ML-206.755 ML-206.756 ML-206.757 ML-206.758 ML-206.759	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.760 ML-206.761 ML-206.762 ML-206.763 ML-206.764 ML-206.765 ML-206.766 ML-206.767	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.768 ML-206.769 ML-206.770 ML-206.771 ML-206.772 ML-206.773 ML-206.774 ML-206.775	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm	



Fig. P3G

ML-206.792 Ø 2 mm X 20 cm ML-206.801 Ø 4 mm X 15 cm ML-206.809 Ø	5 mm X 25 cm 5 mm X 30 cm
ML-206.794       Ø 2 mm X 30 cm       ML-206.803       Ø 4 mm X 25 cm       ML-206.811       Ø 0         ML-206.795       Ø 2 mm X 35 cm       ML-206.804       Ø 4 mm X 30 cm       ML-206.812       Ø 0         ML-206.796       Ø 3 mm X 15 cm       ML-206.805       Ø 4 mm X 35 cm       ML-206.813       Ø 0         ML-206.797       Ø 3 mm X 20 cm       ML-206.806       Ø 5 mm X 15 cm       ML-206.814       Ø 0	5 5 mm X 35 cm 5 6 mm X 35 cm 6 6 mm X 15 cm 6 6 mm X 20 cm 6 6 mm X 25 cm 6 6 mm X 30 cm 6 6 mm X 35 cm








Fig. Py

Power Cannula					
ML-206.831 ML-206.832 ML-206.833 ML-206.834 ML-206.835 ML-206.836 ML-206.837 ML-206.838 ML-206.839	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.840 ML-206.841 ML-206.842 ML-206.843 ML-206.844 ML-206.845 ML-206.847	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.848 ML-206.849 ML-206.850 ML-206.851 ML-206.852 ML-206.853 ML-206.854 ML-206.855	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm



Fig. Le

		Pow	er Cannula		
ML-206.861 ML-206.862 ML-206.863 ML-206.864 ML-206.865 ML-206.866 ML-206.867 ML-206.868 ML-206.869	Ø 2 mm X 15 cm Ø 2 mm X 20 cm Ø 2 mm X 25 cm Ø 2 mm X 30 cm Ø 2 mm X 35 cm Ø 3 mm X 15 cm Ø 3 mm X 20 cm Ø 3 mm X 25 cm Ø 3 mm X 30 cm	ML-206.870 ML-206.871 ML-206.872 ML-206.873 ML-206.874 ML-206.875 ML-206.876 ML-206.877	Ø 3 mm X 35 cm Ø 4 mm X 15 cm Ø 4 mm X 20 cm Ø 4 mm X 25 cm Ø 4 mm X 30 cm Ø 4 mm X 35 cm Ø 5 mm X 15 cm Ø 5 mm X 20 cm	ML-206.878 ML-206.879 ML-206.880 ML-206.881 ML-206.882 ML-206.883 ML-206.884 ML-206.885	Ø 5 mm X 25 cm Ø 5 mm X 30 cm Ø 5 mm X 35 cm Ø 6 mm X 15 cm Ø 6 mm X 20 cm Ø 6 mm X 25 cm Ø 6 mm X 30 cm Ø 6 mm X 35 cm






## Harvester Cannulas



### HARVESTER CANNULAS





**Multiport Speed Harvester- Both End Sharp**A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

	60cc Tommey hub	
ML-301.01	Ø 2.5 mm X 15 cm	
ML-301.02	Ø 2.5 mm X 20 cm	
ML-301.03	Ø 3.0 mm X 15 cm	
ML-301.04	Ø 3.0 mm X 20 cm	
ML-301.05	Ø 3.0 mm X 25 cm	
ML-301.06	Ø 4.0 mm X 25 cm	
ML-301.07	Ø 4.0 mm X 30 cm	
ML-301.08	Ø 4.0 mm X 35 cm	
ML-301.09	Ø 4.5 mm X 30 cm	
ML-301.10	Ø 4.5 mm X 35 cm	

ML-302.01 Ø 3.0 mm X 15 cm ML-302.02 Ø 3.0 mm X 25 cm ML-302.03 Ø 4.0 mm X 15 cm		Treaded Hub	
ML-302.04 Ø 4.0 mm X 25 cm	ML-302.02 ML-302.03	Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm	

	Super Luer Lock	
ML-303.01	Ø 2.0 mm X 15 cm	
ML-303.02	Ø 2.5 mm X 15 cm	
ML-303.03 ML-303.04	Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm	
ML-303.05	Ø 3.0 mm X 20 cm	
ML-303.06	Ø 3.0 mm X 25 cm	


	Luer lock	
ML-304.01	Ø 2.0 mm X 15 cm	
ML-304.02	Ø 2.5 mm X 15 cm	
ML-304.03	Ø 2.5 mm X 20 cm	
ML-304.04	Ø 3.0 mm X 15 cm	
ML-304.05	Ø 3.0 mm X 20 cm	
ML-304.06	Ø 3.0 mm X 25 cm	

Luer lock (Stainless Steel)		
ML-305.01	Ø 2.0 mm X 15 cm	
ML-305.02	Ø 2.5 mm X 15 cm	
ML-305.03	Ø 2.5 mm X 20 cm	
ML-305.04	Ø 3.0 mm X 15 cm	
ML-305.05	Ø 3.0 mm X 20 cm	
ML-305.06	Ø 3.0 mm X 25 cm	

	Power Cannula
ML-306.01	Ø 2 mm X 15 cm
ML-306.02	Ø 2 mm X 20 cm
ML-306.03	Ø 2 mm X 25 cm
ML-306.04	Ø 2 mm X 30 cm
ML-306.05	Ø 2 mm X 35 cm
ML-306.06	Ø 3 mm X 15 cm
ML-306.07	Ø 3 mm X 20 cm
ML-306.08	Ø 3 mm X 25 cm
ML-306.09	Ø 3 mm X 30 cm
ML-306.10	Ø 3 mm X 35 cm
ML-306.11	Ø 4 mm X 15 cm
ML-306.12	Ø 4 mm X 20 cm
ML-306.13	Ø 4 mm X 25 cm
ML-306.14	Ø 4 mm X 30 cm
ML-306.15	Ø 4 mm X 35 cm
ML-306.16	Ø 5 mm X 15 cm
ML-306.17	Ø 5 mm X 20 cm
ML-306.18	Ø 5 mm X 25 cm
ML-306.19	Ø 5 mm X 30 cm
ML-306.20	Ø 5 mm X 35 cm
ML-306.21	Ø 6 mm X 15 cm
ML-306.22	Ø 6 mm X 20 cm
ML-306.23	Ø 6 mm X 25 cm
ML-306.24	Ø 6 mm X 30 cm
ML-306.25	Ø 6 mm X 35 cm



### HARVESTER CANNULAS



**Multiport Speed Harvester- One End Sharp**A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

	60cc Tommey hub	
ML-301.11	Ø 2.5 mm X 15 cm	
ML-301.12	Ø 2.5 mm X 20 cm	
ML-301.13	Ø 3.0 mm X 15 cm	
ML-301.14	Ø 3.0 mm X 20 cm	
ML-301.15	Ø 3.0 mm X 25 cm	
ML-301.16	Ø 4.0 mm X 25 cm	
ML-301.17	Ø 4.0 mm X 30 cm	
ML-301.18	Ø 4.0 mm X 35 cm	
ML-301.19	Ø 4.5 mm X 30 cm	
ML-301.20	Ø 4.5 mm X 35 cm	

	Treaded Hub	
ML-302.11 ML-302.12 ML-302.13 ML-302.14	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	


Super Luer Lock		
ML-303.11	Ø 2.0 mm X 15 cm	
ML-303.12	Ø 2.5 mm X 15 cm	
ML-303.13	Ø 2.5 mm X 20 cm	
ML-303.14	Ø 3.0 mm X 15 cm	
ML-303.15	Ø 3.0 mm X 20 cm	
ML-303.16	Ø 3.0 mm X 25 cm	

	Luer lock
ML-304.11	Ø 2.0 mm X 15 cm
ML-304.12	Ø 2.5 mm X 15 cm
ML-304.13	Ø 2.5 mm X 20 cm
ML-304.14	Ø 3.0 mm X 15 cm
ML-304.15	Ø 3.0 mm X 20 cm
ML-304.16	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)		
ML-305.11 ML-305.12 ML-305.13 ML-305.14 ML-305.15	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-305.16	Ø 3.0 mm X 25 cm	

	Power Cannula
ML-306.11	Ø 2 mm X 15 cm
ML-306.12	Ø 2 mm X 20 cm
ML-306.13	Ø 2 mm X 25 cm
ML-306.14	Ø 2 mm X 30 cm
ML-306.15	Ø 2 mm X 35 cm
ML-306.16	Ø 3 mm X 15 cm
ML-306.17	Ø 3 mm X 20 cm
ML-306.18	Ø 3 mm X 25 cm
ML-306.19	Ø 3 mm X 30 cm
ML-306.20	Ø 3 mm X 35 cm
ML-306.21	Ø 4 mm X 15 cm
ML-306.22	Ø 4 mm X 20 cm
ML-306.23	Ø 4 mm X 25 cm
ML-306.24	Ø 4 mm X 30 cm
ML-306.25	Ø 4 mm X 35 cm
ML-306.26	Ø 5 mm X 15 cm
ML-306.27	Ø 5 mm X 20 cm
ML-306.28	Ø 5 mm X 25 cm
ML-306.29	Ø 5 mm X 30 cm
ML-306.30	Ø 5 mm X 35 cm
ML-306.31	Ø 6 mm X 15 cm
ML-306.32	Ø 6 mm X 20 cm
ML-306.33	Ø 6 mm X 25 cm
ML-306.34	Ø 6 mm X 30 cm
ML-306.35	Ø 6 mm X 35 cm







#### Multi Port Harvester - Both End Bevel

A sophisticated, precision-designed harvester that is truly Patient Friendly. The ports in the Sorensen Harvester are machined with alternating 60 degree cutting edges. This unique design is less traumatic than other multi-port harvesters and delivers a fine-grained, high fat content/purity yield that can easily pass through 1.2mm, 0.9mm and 0.7mm blunt injectors.

	60cc Tommey hub	
ML-301.21	Ø 2.5 mm X 15 cm	
ML-301.22	Ø 2.5 mm X 20 cm	
ML-301.23	Ø 3.0 mm X 15 cm	
ML-301.24	Ø 3.0 mm X 20 cm	
ML-301.25	Ø 3.0 mm X 25 cm	
ML-301.26	Ø 4.0 mm X 25 cm	
ML-301.27	Ø 4.0 mm X 30 cm	
ML-301.28	Ø 4.0 mm X 35 cm	
ML-301.29	Ø 4.5 mm X 30 cm	
ML-301.30	Ø 4.5 mm X 35 cm	

	Treaded Hub	
ML-302.21 ML-302.22 ML-302.23 ML-302.24	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

ML-303.21 Ø 2.0 mm X 15 cm ML-303.22 Ø 2.5 mm X 15 cm ML-303.23 Ø 2.5 mm X 20 cm ML-303.24 Ø 3.0 mm X 15 cm ML-303.25 Ø 3.0 mm X 20 cm		Super Luer Lock	
WIL-303.20 Ø 3.0 HIIII X 23 GH	ML-303.22 ML-303.23 ML-303.24	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm	

	Luer lock	
ML-304.21	Ø 2.0 mm X 15 cm	
ML-304.22 ML-304.23	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm	
ML-304.24	Ø 3.0 mm X 15 cm	
ML-304.25	Ø 3.0 mm X 20 cm	
ML-304.26	Ø 3.0 mm X 25 cm	

	Luer lock (Stainless Steel)
ML-305.21 ML-305.22 ML-305.23 ML-305.24 ML-305.25	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm
ML-305.26	Ø 3.0 mm X 25 cm

	Power Cannula
ML-306.21	Ø 2 mm X 15 cm
ML-306.22	Ø 2 mm X 20 cm
ML-306.23	Ø 2 mm X 25 cm
ML-306.24	Ø 2 mm X 30 cm
ML-306.25	Ø 2 mm X 35 cm
ML-306.26	Ø 3 mm X 15 cm
ML-306.27	Ø 3 mm X 20 cm
ML-306.28	Ø 3 mm X 25 cm
ML-306.29	Ø 3 mm X 30 cm
ML-306.30	Ø 3 mm X 35 cm
ML-306.31	Ø 4 mm X 15 cm
ML-306.32	Ø 4 mm X 20 cm
ML-306.33	Ø 4 mm X 25 cm
ML-306.34	Ø 4 mm X 30 cm
ML-306.35	Ø 4 mm X 35 cm
ML-306.36	Ø 5 mm X 15 cm
ML-306.37	Ø 5 mm X 20 cm
ML-306.38	Ø 5 mm X 25 cm
ML-306.39	Ø 5 mm X 30 cm
ML-306.40	Ø 5 mm X 35 cm
ML-306.41	Ø 6 mm X 15 cm
ML-306.42	Ø 6 mm X 20 cm
ML-306.43	Ø 6 mm X 25 cm
ML-306.44	Ø 6 mm X 30 cm
ML-306.45	Ø 6 mm X 35 cm





## Injector Cannulas

(Micro Injectors)



### MICRO INJECTORS















This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

	Super Luer Lock		Luer lock
ML-403.02 ML-403.03 ML-403.04	Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 3.0 mm X 15 cm	ML-404.01 ML-404.02 ML-404.03 ML-404.04 ML-404.05	Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)						
3cm 5cm 7cm 10cm 12cm						
1.1Ø 1.4Ø 1.8Ø 2.0Ø	ML-405.01 ML-405.02 ML-405.03 ML-405.04	ML-405.11 ML-405.12 ML-405.13 ML-405.14	ML-405.21 ML-405.22 ML-405.23 ML-405.24	ML-405.31 ML-405.32 ML-405.33 ML-405.34	ML-405.41 ML-405.42 ML-405.43 ML-405.44	

#### Blunt Extractor/Injector

This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

	Super Luer Lock		Luer lock
ML-403.11	Ø 2.0 mm X 15 cm	ML-404.11	Ø 2.0 mm X 15 cm
ML-403.12	Ø 2.5 mm X 15 cm	ML-404.12	Ø 2.5 mm X 15 cm
ML-403.13	Ø 3.0 mm X 15 cm	ML-404.13	Ø 3.0 mm X 15 cm

	Luer lock (Stainless Steel)						
3cm 5cm 7cm 10cm 12cm							
	1.1Ø 1.4Ø 1.8Ø 2.0Ø	ML-405.51 ML-405.52 ML-405.53 ML-405.54	ML-405.61 ML-405.62 ML-405.63 ML-405.64	ML-405.71 ML-405.72 ML-405.73 ML-405.74	ML-405.81 ML-405.82 ML-405.83 ML-405.84	ML-405.91 ML-405.92 ML-405.93 ML-405.94	

**Blunt Tip Injector**Rounded blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp instruments.

	Super Luer Lock		Luer lock
ML-403.21	Ø 1.0 mm X 4 cm	ML-404.21	Ø 1.0 mm X 4 cm
ML-403.22	Ø 1.0 mm X 5 cm	ML-404.22	Ø 1.0 mm X 5 cm
ML-403.23	Ø 1.0 mm X 7 cm	ML-404.23	Ø 1.0 mm X 7 cm
ML-403.24	Curved 12 cm	ML-404.24	Curved 12 cm
ML-403.25	Curved 15 cm	ML-404.25	Curved 15 cm
ML-403.26	Curved 20 cm	ML-404.26	Curved 20 cm

Luer lock (Stainless Steel)							
3cm 5cm 7cm 10cm 12cm							
1.1Ø 1.4Ø 1.8Ø	ML-405.102	ML-405.111 ML-405.112 ML-405.113	ML-405.122	ML-405.132	ML-405.142		
	ML-405.103						

## MICRO INJECTORS





	Luer lock (Stainless Steel)							
		3cm	5cm	7cm	10cm	12cm		
1	l.1Ø	ML-405.201	ML-405.211	ML-405.221	ML-405.231	ML-405.241		
1	1.4Ø	ML-405.202	ML-405.212	ML-405.222	ML-405.232	ML-405.242		
1	l.8Ø	ML-405.203	ML-405.213	ML-405.223	ML-405.233	ML-405.243		
2	2.0Ø	ML-405.204	ML-405.214	ML-405.224	ML-405.234	ML-405.244		

Fig. 1C

	Luer lock (Stainless Steel)						
	3cm	5cm	7cm	10cm	12cm		
1.1Ø			ML-405.271				
1.4Ø		ML-405.262					
1.8Ø		ML-405.263					
2.0Ø	ML-405.254	ML-405.264	ML-405.274	ML-405.284	ML-405.294		


Fig. 1P

	Luer lock (Stainless Steel)							
	3cm	5cm	7cm	10cm	12cm			
1.1Ø 1.4Ø 1.8Ø 2.0Ø	ML-405.302 ML-405.303	ML-405.312	ML-405.323	ML-405.332 ML-405.333	ML-405.342 ML-405.343			

	Luer lock (Stainless Steel)							
3cm			5cm	7cm	10cm	12cm		
	1.1Ø		ML-405.361			ML-405.391		
	1.4Ø	ML-405.352	ML-405.362	ML-405.372	ML-405.382	ML-405.392		
	1.8Ø	ML-405.353	ML-405.363	ML-405.373	ML-405.383	ML-405.393		
	2.0Ø	ML-405.354	ML-405.364	ML-405.374	ML-405.384	ML-405.394		







### Coleman Liposuction Cannula Fig. MIC I

	Luer lock (Stainless Steel)									
		3cm	5cm	7cm	10cm	12cm	15cm	20cm	3.5cm	
-	1.1Ø	ML-405.401	ML-405.411	ML-405.421	ML-405.431	ML-405.441	ML-405.451	ML-405.461	ML-405.471	
	1.4Ø	ML-405.402	ML-405.412	ML-405.422	ML-405.432	ML-405.442	ML-405.452	ML-405.462	ML-405.472	
	1.8Ø	ML-405.403	ML-405.413	ML-405.423	ML-405.433	ML-405.443	ML-405.453	ML-405.463	ML-405.473	
2	2.0Ø	ML-405.404	ML-405.414	ML-405.424	ML-405.434	ML-405.444	ML-405.454	ML-405.464	ML-405.474	
2	2.5Ø	ML-405.405	ML-405.415	ML-405.425	ML-405.435	ML-405.445	ML-405.455	ML-405.465	ML-405.475	

### Coleman Liposuction Cannula Fig. MIC II

		Luer lock (Stainless Steel)								
		3cm	5cm	7cm	10cm	12cm	15cm	20cm		
ı	1.1Ø	ML-405.501	ML-405.511	ML-405.521	ML-405.531	ML-405.541	ML-405.551	ML-405.561		
	1.4Ø	ML-405.502	ML-405.512	ML-405.522	ML-405.532	ML-405.542	ML-405.552	ML-405.562		
	1.8Ø	ML-405.503	ML-405.513	ML-405.523	ML-405.533	ML-405.543	ML-405.553	ML-405.563		
	2.0Ø	ML-405.504	ML-405.514	ML-405.524	ML-405.534	ML-405.544	ML-405.554	ML-405.564		
	2.5Ø	ML-405.505	ML-405.515	ML-405.525	ML-405.535	ML-405.545	ML-405.555	ML-405.565		



## **Coleman Liposuction Cannula Spoon Tip Fig. MIC III**This tip works well for needle extraction of soft tissue and as a soft tissue re injection needle.

	Luer lock (Stainless Steel)							
	3cm	5cm	7cm	10cm	12cm	15cm	20cm	
1.1Ø	ML-405.601	ML-405.611	ML-405.621	ML-405.631	ML-405.641	ML-405.651	ML-405.661	
1.4Ø	ML-405.602	ML-405.612	ML-405.622	ML-405.632	ML-405.642	ML-405.652	ML-405.662	
1.8Ø	ML-405.603	ML-405.613	ML-405.623	ML-405.633	ML-405.643	ML-405.653	ML-405.663	
2.0Ø	ML-405.604	ML-405.614	ML-405.624	ML-405.634	ML-405.644	ML-405.654	ML-405.664	
2.5Ø	ML-405.605	ML-405.615	ML-405.625	ML-405.635	ML-405.645	ML-405.655	ML-405.665	



#### **Little Injector**

Open and more pointed blunt tip is very effective in dense and fibrous tissue.

	Super Luer Lock	
ML-403.51	Ø 1.0 mm X 4 cm	ML
ML-403.52 ML-403.53	Ø 1.0 mm X 5 cm Ø 1.0 mm X 7 cm	ML ML

Luer lock				
ML-404.51	Ø 1.0 mm X 4 cm			
ML-404.52	Ø 1.0 mm X 5 cm			
ML <b>-</b> 404.53	Ø 1.0 mm X 7 cm			

	Luer lock (Stainless Steel)								
	3cm	5cm	7cm	10cm	12cm				
1.10 1.40 1.80 2.00	ML-405.702 ML-405.703	ML-405.712 ML-405.713	ML-405.721 ML-405.722 ML-405.723 ML-405.724	ML-405.732 ML-405.733	ML-405.741 ML-405.742 ML-405.743 ML-405.744				

### MICRO INJECTORS







Fig. 1B

	Luer lock (Stainless Steel)						
	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	ML-405.801	ML-405.811	ML-405.821	ML-405.831	ML-405.841	ML-405.851	ML-405.861
1.4Ø	ML-405.802	ML-405.812	ML-405.822	ML-405.832	ML-405.842	ML-405.852	ML-405.862
1.8Ø	ML-405.803	ML-405.813	ML-405.823	ML-405.833	ML-405.843	ML-405.853	ML-405.863
2.0Ø	ML-405.804	ML-405.814	ML-405.824	ML-405.834	ML-405.844	ML-405.854	ML-405.864
2.5Ø	ML-405.805	ML-405.815	ML-405.825	ML-405.835	ML-405.845	ML-405.855	ML-405.865

### **Curved Injectors**



#### **Curved Injectors**

New Curved Injectors are designed especially for fat transfer procedures to the body. Their unique shape mimics the anatomical curvature of the female torso and, thus, easily produces the natural-looking results that patients desire. Curved Injection Cannulas ameliorate the fat transfer process by providing a direct and precise route to the injection site and by lessening the time required to correctly sculpt curves.






# Speciality Cannulas

Speciality Cannula are useful tools for different to operation are such as face or Neck.



### SPECIALITY CANNULAS





#### **V- Dissector**

This cannula is specifically designed for use in the superficial compartment of the truncal area. The cutting notch is very useful for undermining in facial procedures and in dense, fibrous tissue.

	60cc Tommey hub	
ML-501.01	Ø 2.5 mm X 15 cm	
ML-501.02	Ø 2.5 mm X 20 cm	
ML-501.03	Ø 3.0 mm X 15 cm	
ML-501.04	Ø 3.0 mm X 20 cm	
ML-501.05	Ø 3.0 mm X 25 cm	
ML-501.06	Ø 4.0 mm X 25 cm	
ML-501.07	Ø 4.0 mm X 30 cm	
ML-501.08	Ø 4.0 mm X 35 cm	
ML-501.09	Ø 4.5 mm X 30 cm	
ML-501.10	Ø 4.5 mm X 35 cm	

	Treaded Hub	
ML-502.01 ML-502.02 ML-502.03 ML-502.04	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

	Super Luer Lock	
ML-503.01 ML-503.02 ML-503.03 ML-503.04 ML-503.05 ML-503.06	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm	


	Luer lock	
ML-504.01	Ø 2.0 mm X 15 cm	
ML-504.02	Ø 2.5 mm X 15 cm	
ML-504.03	Ø 2.5 mm X 20 cm	
ML-504.04	Ø 3.0 mm X 15 cm	
ML-504.05	Ø 3.0 mm X 20 cm	
ML-504.06	Ø 3.0 mm X 25 cm	

Luer lock (Stainless Steel)		
ML-505.01	Ø 2.0 mm X 15 cm	
ML-505.02 ML-505.03	Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm	
ML-505.04 ML-505.05	Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm	
ML-505.06	Ø 3.0 mm X 25 cm	

	Power Cannula
ML-506.01	Ø 2 mm X 15 cm
ML-506.02	Ø 2 mm X 20 cm
ML-506.03	Ø 2 mm X 25 cm
ML-506.04	Ø 2 mm X 30 cm
ML-506.05	Ø 2 mm X 35 cm
ML-506.06	Ø 3 mm X 15 cm
ML-506.07	Ø 3 mm X 20 cm
ML-506.08	Ø 3 mm X 25 cm
ML-506.09	Ø 3 mm X 30 cm
ML-506.10	Ø 3 mm X 35 cm
ML-506.11	Ø 4 mm X 15 cm
ML-506.12	Ø 4 mm X 20 cm
ML-506.13	Ø 4 mm X 25 cm
ML-506.14	Ø 4 mm X 30 cm
ML-506.15	Ø 4 mm X 35 cm
ML-506.16	Ø 5 mm X 15 cm
ML-506.17	Ø 5 mm X 20 cm
ML-506.18	Ø 5 mm X 25 cm
ML-506.19	Ø 5 mm X 30 cm
ML-506.20	Ø 5 mm X 35 cm
ML-506.21	Ø 6 mm X 15 cm
ML-506.22	Ø 6 mm X 20 cm
ML-506.23	Ø 6 mm X 25 cm
ML-506.24	Ø 6 mm X 30 cm
ML-506.25	Ø 6 mm X 35 cm









#### Flap

The Flap is primarily used as an undermining dissector, and in smaller diameters, for facial tissue removal. The Flap Cannula in a 2.1 mm diameter is very useful for undermining the nasolabial fold prior to autologous tissue injections, as well as for tissue removal in this anatomic area.

injections, as w	vell as for tissue removal in this anatomic
	60cc Tommey hub
ML-501.11	Ø 2.5 mm X 15 cm
ML-501.12	Ø 2.5 mm X 20 cm
ML-501.13 ML-501.14	Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm
ML-501.14 ML-501.15	Ø 3.0 mm X 25 cm
ML-501.16	Ø 4.0 mm X 25 cm
ML-501.17	Ø 4.0 mm X 30 cm
ML-501.18	Ø 4.0 mm X 35 cm
ML-501.19	Ø 4.5 mm X 30 cm
ML-501.20	Ø 4.5 mm X 35 cm
	Treaded Hub
	Treatieu nub
ML-502.11	Ø 3.0 mm X 15 cm
ML-502.12 ML-502.13	Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm
ML-502.13	Ø 4.0 mm X 25 cm
WE 002.11	2 110 11111 7 20 0111
	Super Luer Lock
ML-503.11	Ø 2.0 mm X 15 cm
ML-503.12	Ø 2.5 mm X 15 cm
ML-503.13	Ø 2.5 mm X 20 cm
ML-503.14	Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm
ML-503.15 ML-503.16	Ø 3.0 mm X 20 cm
WIE-503. 10	D J.O HIM A ZJ GH
	Luer lock

	Luer lock	
ML-504.11	Ø 2.0 mm X 15 cm	
ML-504.12	Ø 2.5 mm X 15 cm	
ML-504.13	Ø 2.5 mm X 20 cm	
ML-504.14	Ø 3.0 mm X 15 cm	
ML-504.15	Ø 3.0 mm X 20 cm	
ML-504.16	Ø 3.0 mm X 25 cm	

	Luer lock (Stainless Steel)
ML-505.11 ML-505.12 ML-505.13 ML-505.14 ML-505.15 ML-505.16	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm Ø 2.5 mm X 20 cm Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm Ø 3.0 mm X 25 cm
WE 000.10	Power Connuls

	Power Cannula
ML-506.51	Ø 2 mm X 15 cm
ML-506.52	Ø 2 mm X 20 cm
ML-506.53	Ø 2 mm X 25 cm
ML-506.54	Ø 2 mm X 30 cm
ML-506.55	Ø 2 mm X 35 cm
ML-506.56	Ø 3 mm X 15 cm
ML-506.57	Ø 3 mm X 20 cm
ML-506.58	Ø 3 mm X 25 cm
ML-506.59	Ø 3 mm X 30 cm
ML-506.60	Ø 3 mm X 35 cm
ML-506.61	Ø 4 mm X 15 cm
ML-506.62	Ø 4 mm X 20 cm
ML-506.63	Ø 4 mm X 25 cm
ML-506.64	Ø 4 mm X 30 cm
ML-506.65	Ø 4 mm X 35 cm
ML-506.66	Ø 5 mm X 15 cm
ML-506.67	Ø 5 mm X 20 cm
ML-506.68	Ø 5 mm X 25 cm
ML-506.69	Ø 5 mm X 30 cm
ML-506.70	Ø 5 mm X 35 cm
ML-506.71	Ø 6 mm X 15 cm
ML-506.72	Ø 6 mm X 20 cm
ML-506.73	Ø 6 mm X 25 cm
ML-506.74	Ø 6 mm X 30 cm
ML-506.75	Ø 6 mm X 35 cm



### SPECIALITY CANNULAS



### Field Flap

With its elevated leading edge, this cannula is designed to open and expose soft tissue for more precise suctioning along the mandible and sub-mental areas.

	60cc Tommey hub
ML-501.21	Ø 2.5 mm X 15 cm
ML-501.22	Ø 2.5 mm X 20 cm
ML-501.23	Ø 3.0 mm X 15 cm
ML-501.24	Ø 3.0 mm X 20 cm
ML-501.25	Ø 3.0 mm X 25 cm
ML-501.26	Ø 4.0 mm X 25 cm
ML-501.27	Ø 4.0 mm X 30 cm
ML-501.28	Ø 4.0 mm X 35 cm
ML-501.29	Ø 4.5 mm X 30 cm
ML-501.30	Ø 4.5 mm X 35 cm

	Treaded Hub	
ML-502.21 ML-502.22 ML-502.23 ML-502.24	Ø 3.0 mm X 15 cm Ø 3.0 mm X 25 cm Ø 4.0 mm X 15 cm Ø 4.0 mm X 25 cm	

	Super Luer Lock	
ML-503.21	Ø 2.0 mm X 15 cm	
ML-503.22	Ø 2.5 mm X 15 cm	
ML-503.23	Ø 2.5 mm X 20 cm	
ML-503.24	Ø 3.0 mm X 15 cm	
ML-503.25	Ø 3.0 mm X 20 cm	
ML-503.26	Ø 3.0 mm X 25 cm	

	Luer lock
ML-504.21 ML-504.22	Ø 2.0 mm X 15 cm Ø 2.5 mm X 15 cm
ML-504.23	Ø 2.5 mm X 20 cm
ML-504.24 ML-504.25	Ø 3.0 mm X 15 cm Ø 3.0 mm X 20 cm
ML-504.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)				
ML-505.21	Ø 2.0 mm X 15 cm			
ML-505.22	Ø 2.5 mm X 15 cm			
ML-505.23	Ø 2.5 mm X 20 cm			
ML-505.24	Ø 3.0 mm X 15 cm			
ML-505.25	Ø 3.0 mm X 20 cm			
ML-505.26	Ø 3.0 mm X 25 cm			

Power Cannula				
ML-506.101	Ø 2 mm X 15 cm			
ML-506.102	Ø 2 mm X 20 cm			
ML-506.103	Ø 2 mm X 25 cm			
ML-506.104	Ø 2 mm X 30 cm			
ML-506.105	Ø 2 mm X 35 cm			
ML-506.106	Ø 3 mm X 15 cm			
ML-506.107	Ø 3 mm X 20 cm			
ML-506.108	Ø 3 mm X 25 cm			
ML-506.109	Ø 3 mm X 30 cm			
ML-506.110	Ø 3 mm X 35 cm			
ML-506.111	Ø 4 mm X 15 cm			
ML-506.112	Ø 4 mm X 20 cm			
ML-506.113	Ø 4 mm X 25 cm			
ML-506.114	Ø 4 mm X 30 cm			
ML-506.115	Ø 4 mm X 35 cm			
ML-506.116	Ø 5 mm X 15 cm			
ML-506.117	Ø 5 mm X 20 cm			
ML-506.118	Ø 5 mm X 25 cm			
ML-506.119	Ø 5 mm X 30 cm			
ML-506.120	Ø 5 mm X 35 cm			
ML-506.121	Ø 6 mm X 15 cm			
ML-506.122	Ø 6 mm X 20 cm			
ML-506.123	Ø 6 mm X 25 cm			
ML-506.124	Ø 6 mm X 30 cm			
ML-506.125	Ø 6 mm X 35 cm			

### SINGLE-USE INSTRUMENTS





#### Super Luer Lock

The Super Luer Lock is on our Single-use instruments.

#### **Tumescent Infiltrator**

This cannula is used specifically for the infusion of tumescent solution and for hydraulic dissection. the configuration of the ports allows even distribution of the fluid.

SL-601.01 Ø 2.0 mm X 12cm SL-601.02 Ø 2.0 mm X 20cm



#### Tonnard Harvester

A multiport harvester with 1mm elevated sharp microports. The 20 port configuration allows for harvesting both on the forward and backward strokes.

SL-602.01 Ø 3.0 mm X 20cm SL-602.02 Ø 2.4 mm X 20cm



#### Tonnard Harvester

Excellent for harvesting small aliquots of tissues.

SL-603.01 Ø 2.5 mm X 15cm SL-603.02 Ø 2.0 mm X 15cm



#### Carraway Harvester

Spiral three port design effective for multi-plane tissue harvesting.

SL-604.01 Ø 2.5 mm X 15cm SL-604.02 Ø 2.0 mm X 10cm SL-604.03 Ø 2.0 mm X 15cm SL-604.04 Ø 1.5 mm X 8cm


SL-604.05 Ø 1.5 mm X 10cm



### Injector

Round blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp injectors.

SL-605.01 Ø 1.5 mm X 9cm SL-605.02 Ø 1.5 mm X 7cm SL-605.03 Ø 1.0 mm X 5cm SL-605.04 Ø 1.0 mm X 4cm





### Spoon Tip Injector

Blunt injector designed with a more open port for rapid injection of larger volumes - a torso and buttock workhourse.

SL-606.01 Ø 3.0 mm X 20cm SL-606.02 Ø 2.0 mm X 15cm



#### Miller Speed Harvester

SL-607.01 Ø 2.0 mm X 15cm



## Liposuction Accessories

### Handles

### Connectors

- a) Luer to Luer
- b) Luer to Handle

### **Transfer Adopters**

- a) Luer to Luer 1.2 mm Luer to Luer 1.4 mm Luer to Luer 2.4 mm
- b) Anaerobic Transfer
- c) Micro to Nano

### Snap Lock

- a) Standard Snap Locks 60cc, 20cc, 10cc
- b) Snap Lock

Fat Transfer Gun

**Decanting Stands** 

Sterilization Boxed



Anaerobic Transfer

### **HANDLES FOR CANNULAS**







### HANDLES FOR CANNULAS

ML-704.01 with connection thread for cannulas hose connection



ML-705.01

with connection thread for cannulas hose connection



with connection thread for cannulas hose connection



ML-707.01

with connection thread for cannulas hose connection



ML-750.01

Trocar



The Pocar is a percutaneous punch that is used to make a eussit dna noitartlfini ot roirp noisicni laitnerefmucric extraction. It is designed to maintain a tight seal while minimizing scarring.



### Anaerobic Transfer

### NanoTransfer

### ML-751.01

NanoTransfer



#### ML-752.01

Mesh for NanoTransfer





#### **Anaerobic Transfers**

ML-760.01 Luer 1.2mm ML-760.02 Luer 1.4mm ML-760.03 Luer 2.4mm







#### ML-771.01

Luer to Luer



#### **Decant Stands**

ML-772.02 60cc to Luer ML-773.03 Luer to Luer



### SnapLoks

These syringe snap locks fit into the lunger flanges of the syringe. When the plunger is pulled out the lock snaps open on the barrel lip thus holding the syringe under vacuum.

ML-901.01	60cc	
ML-901.02	50cc	
ML-901.03	30cc	
ML-901.04	20cc	
ML-901.05	10cc	







**Universal Power Injector** 

ML-801.01 Universal Fat Injector 10cc & 20cc Two in One

ML-801.02 Fat Injector 60cc



### **Fat Syringe**

ML-950.01

The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or a 60cc Monoject syringe.







### **Syringe Decanting and Injection**

60 cc
50 cc
30 cc
20 cc
10 cc

medastepro@gmail.com





# Liposuction

Cannulas & Accessories

